

Disposable gown is made from high-performance isolation protective material.

Widely used in medical safety protective, aseptic workshop, protective isolation, mining electronics factory, food factory farm animal husbandry biohazard and so on.

With two kinds processes: sewing and ultrasound thermal bonding.

Material features :

1. PP disposable gown

It is made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Color: white, green, blue, yellow, orange and so on

Material weight: 20-65gsm.

2. PP+PE disposable gown

It is made from PP+PE material, Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

Color: white, green, blue, yellow.

Material weight: 40-65gsm.

3. SMS disposable gown (Standard and Reinforced):

It is made from hydrophobic SMS/PP/PP+PE material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: white, blue.

Material weight: 35-65gsm.

Disposable gown style:

Tie/ hook and loop fastener on collar, tie on waist, elastic/knitted cuff

Disposable gown Specifications

|             |   |
|-------------|---|
| Description | disposable level1 level2 level3 isolation gowns SMS/PP+PE/PP protective cloth |
| Material    | PP, PP+PE, SMS, CPE   |

|                |  |
|----------------|--|
| Weight         | PP 20-65gsm; PP+PE 40-65 gsm; SMS 35-65 gsm (comfortable,without glass fibres,latex free.)   |
| Style          | tie/ hook and loop fastener on collar,tie on waist,elastic/knitted cuff  |
| Size           | 120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm   |
| Color          | white,green,blue,yellow  |
| packing        | non sterile,10 pcs/bag,10 bags/ctn.<br>sterile,1 pcs/sterile bag,100 pcs/ctn   |
| packing design | all inner box and carton printing can be as your design  |
| carton size    | 58X28X48 cm or as your requirement   |
| Applications   | hospital clinical medical personnel and patients.<br>dust-free workshop,laboratory, food industry, electronic manufactures and so on |
| MOQ            | 10000 pcs  |
| Sample         | can be supplied free for your quality checking within 3 days   |
| Payment term   | T/T,PayPal   |
| Delivery time  | as your quantity(we can finish 20,000 pcs per week)  |
| FOB port       | Wuhan/Shanghai/Ningbo/Guangzhou  |



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

