

la especificación de Sentimiento fresco Máscara facial desechable

Tela MATERIAL

Capa 1: spunlace súper suave

Capa 2: Papel de alta filtración

Layer3: Material ES

Estilo

Bucle oído

Dimensiones S IZE

17.5 * 9,5 cm

Tela W OCHO

25GSM + 25GSM + 35GSM, C Ould SER Como su solicitud

color

W HITE, bebé azul, bebé PINK, gris

Certificado

CE, ISO13485, FDA

MOQ:

100000pcs

Diseño de embalaje

Puede ser como su diseño

Embalaje

50 PCS / B OX, 200 0 PC S / CTN, OR Como su requisito

Capacidad

60,000 PCS por día

El tiempo de entrega

7-30 días


APLICACIONES

HOSPITAL, Laboratorio, clínica, Industria alimentaria, Taller sin polvo


la imagen de Sentimiento fresco Máscara facial desechable



Precision welding
resistant to pull-off



Antiparticulate
microorganism




Anti-virus, anti-bacteria



orient Honest
MEDICAL INTERNATIONAL

0.32 cool sense of value



test item			
instant cool value	0.2-0.3	0.3-0.4	>0.4
somatosensory	not obvious	cool without ice	Have a cold feeling

* Data from the laboratory, different colors, fabric batches will have a certain impact on the value