

La propriété de Disposable White PP lab coat

Tout cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: S-5XL or customized
2. Made of nonwoven Matériel, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

The Specification of Disposable White PP lab coat

Nom	Disposable White PP lab coat
Matériel	pp
Style	Tricoté collar and knit cuffs, with bottons and pockets
Taille	S-5XL
Couleur	bleu
emballage	5pcs/box,10 bags/ctn
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Échantillon	can be supplied free for your quality checking within 2 jours
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Délai de livraison	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 jours;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Disposable White PP lab coat





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

