

La spécification de Sentiment cool masque facial jetable

en tissu Matériel

Couche 1: Spunlace super douce

Couche 2: papier de filtration élevé

Layer3: Matériau es

Style

Boucle d'oreille

Dimensions s ize

17.5 * 9,5 cm

Tissu w huit

25gsm + 25gsm + 35gsm, C Ould BE Comme votre demande

couleur

w Hite, de bébé bleu, de bébé rose, gris

Certificat

CE, ISO13485, FDA

MOQ:

100000pcs

Conception d'emballage

peut être comme votre conception

Emballage

50 PCS / B OX, 200 0 PC S / CTN, ou comme votre exigence

Aptitude

60 000 PC par jour

Heure de livraison

7-30 jours


Applications

Hôpital, Laboratoire, clinique, Industrie alimentaire, Atelier sans poussière

L'image de **Sentiment cool** masque facial jetable



Precision welding
resistant to pull-off



Antiparticulate
microorganism





Anti-virus, anti-bacteria



orient Honest
MEDICAL INTERNATIONAL

0.32 cool sense of value



test item			
instant cool value	0.2-0.3	0.3-0.4	>0.4
somatosensory	not abvious	cool eithout ice	Have a cold feeling

* Data from the laboratory, different colors, fabric batches will have a certain impact on the value