

La spécification de Sentiment cool masque facial jetable

en tissu Matériel

Couche 1: Spunlace super douce

Couche 2: papier de filtration élevé

Layer3: Matériau es

Style

Boucle d'oreille

Dimensions s ize

17.5 * 9,5 cm

Tissu w huit

25gsm + 25gsm + 35gsm, C Ould BE Comme votre demande

couleur

w Hite, de bébé bleu, de bébé rose, gris

Certificat

CE, ISO13485, FDA

MOQ:

100000pcs

Conception d'emballage

peut être comme votre conception

Emballage

50 PCS / B OX, 200 0 PC S / CTN, ou comme votre exigence

Aptitude

60 000 PC par jour

Heure de livraison

7-30 jours

Applications

Hôpital, Laboratoire, clinique, Industrie alimentaire, Atelier sans poussière

L'image de Sentiment cool masque facial jetable




One more layer of protection,
one more peace of mind





orient Honest
MEDICAL INTERNATIONAL

0.32 cool sense of value



test item			
instant cool value	0.2-0.3	0.3-0.4	>0.4
somatosensory	not obvious	cool without ice	Have a cold feeling

* Data from the laboratory, different colors, fabric batches will have a certain impact on the value