

## **La spécification de Sentiment cool masque facial jetable**

### **en tissu Matériel**

Couche 1: Spunlace super douce

Couche 2: papier de filtration élevé

Layer3: Matériau es

### **Style**

Boucle d'oreille

### **Dimensions s ize**

17.5 \* 9,5 cm

### **Tissu w huit**

25gsm + 25gsm + 35gsm, C Ould BE Comme votre demande

### **couleur**

w Hite, de bébé bleu, de bébé rose, gris

### **Certificat**

CE, ISO13485, FDA

### **MOQ:**

100000pcs

### **Conception d'emballage**

peut être comme votre conception

### **Emballage**

50 PCS / B OX, 200 0 PC S / CTN, ou comme votre exigence

### **Aptitude**

60 000 PC par jour

### **Heure de livraison**

7-30 jours

### **Applications**

Hôpital, Laboratoire, clinique, Industrie alimentaire, Atelier sans poussière

## L'image de Sentiment cool masque facial jetable




One more layer of protection,  
one more peace of mind





orient Honest  
MEDICAL INTERNATIONAL

## 0.32 cool sense of value



test item			
instant cool value	0.2-0.3	0.3-0.4	>0.4
somatosensory	not obvious	cool without ice	Have a cold feeling

\* Data from the laboratory, different colors, fabric batches will have a certain impact on the value