

**Berkualitas tinggi sekali pakai non anyaman penutup kasur**

Gaya: datar, ujung elastis, ujung persegi

Ukuran: 80x190+20cm, 90x200+20cm, 160x190+20cm, 160x200+20cm, 180x200+20cm ...

Pengemasan: 5 pcs/tas, 5pcs/ag atau 1 pcs/tas, 50 pcs/ctn.

Bahan kain	Pp, pp/pe, sms, smms,
Gaya	Datar, ujung elastis, ujung persegi, gulung
Ukuran dimensi	80x190+20cm, 90x200+20cm, 160x190+20cm, 160x200+20cm, 180x200+20cm
Warna	putih, hijau, biru, kuning, merah muda dll
Sertifikat	CE, ISO13485, FDA
MOQ:	5000 PC
Desain Pengepakan	Semua tas dalam dan pencetakan karton bisa menjadi desain Anda
Sedang mengemas	10 pcs/box, 10 tas/ctn atau 5 pcs/tas, 50 pcs/tas
Kemampuan	Kapasitas: 10.000 pcs/hari,
Waktu pengiriman	1*40HQ Lembar tidur sekali pakai dapat diselesaikan dalam waktu 10-15 hari;
Aplikasi	Rumah Sakit, Laboratorium, Klinik, Industri Makanan, Lokakarya Bebas Debu
Port fob	Wuhan/Shanghai/Ningbo/Guangzhou
Cara pembayaran	T/t, paypal, western union, l/c, d/a, d/p, escrow,



orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL DISTRIBUTOR





orient Honest  
MEDICAL TEXTILES

