

**All isolation gowns are made of the highest quality materials, providing excellent protection and superior comfort:**

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Fabric Material</b>     | PP+PE material                                                 |
| <b>Style</b>               | With ties on collar and waist, knitted cuffs                   |
| <b>Dimensions size</b>     | X-3XL                                                          |
| <b>Fabric weight</b>       | 30-40g, could be as your request                               |
| <b>Color</b>               | Green color                                                    |
| <b>Certificate</b>         | ISO13485, FDA                                                  |
| <b>MOQ:</b>                | 10000 pcs                                                      |
| <b>Packing design</b>      | can be as your design                                          |
| <b>Packing</b>             | 1pc/plastic paper bag,50pcs/ctn, or as your requirement        |
| <b>Capability</b>          | 60,000pcs per day                                              |
| <b>Delivery time</b>       | 7-30 days                                                      |
| <b><u>Applications</u></b> | laboratory, clinic, food industry, dust-free workshop,hospital |

**The picture of disposable Level 1 isolation gown with knitted cuffs medical manufacturer**







