

### Properti dari Disposable SMS lab coat

**Semua cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:**

- 1.Size: S-5XL or customized
2. Made of nonwoven bahan, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

### The Specification of Sekali pakai SMS Jas laboratorium

Nama	Sekali pakai SMS Jas laboratorium
Bahan	SMS
Gaya	Knitted collar and knit cuffs, with bottons and pockets
Ukuran	S-5XL
Warna	biru
sedang mengemas	5pcs/box,10 bags/ctn
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Sampel	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

### Properti dari Disposable SMS lab coat

**Semua cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:**

- 1.Size: S-5XL or customized
2. Made of nonwoven bahan, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

### The Specification of Sekali pakai SMS Jas laboratorium

Nama	Sekali pakai SMS Jas laboratorium
Bahan	SMS
Gaya	Knitted collar and knit cuffs, with bottons and pockets
Ukuran	S-5XL
Warna	biru
sedang mengemas	5pcs/box,10 bags/ctn
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Sampel	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL







orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

