

### Key Specifications/Special Features:

Material: PP nonwoven

Gram: 40gsm or customized

Size: general size or customized

Color: blue/pink/black, etc

Application: SPA, Massage, Hospital, Salon, etc

Storage: stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

Sample: free

Printed Logo: Yes

### Feature:

1. Strict standard quality control with CE, ISO.
2. Made of high quality of PP material.
3. Factory direct sales with competitive price.
4. Experienced staff, fast delivery, stable production capacity.
5. Over 10 years production experience.
6. OEM is available, different sizes, thickness, colors, printed logos, etc.
7. Professional sales teams to serve you wholeheartedly.

### Packing & Shipping:

1. Normal Package: 1pcs/bag, 5pcs/mid bag, 600pcs/carton.
2. MOQ: 10000pcs.
3. Port: Shanghai / Wuhan, China.
4. Sample availability: YES.
5. Normal lead time: 3 days for sample order, 25 days for bulk order.

### Why Choose Us:

1. We have passed CE, FDA, ISO13485 certification.
2. Welcome any OEM requirements.
3. Advanced equipment and professional staffs with over 10 years export experience.
4. Qualified products, 100% new brand material, safe and sanitary.
5. We have strict QC quality control system.
6. Have enough experience and gained a high reputation from our customers all the world. Such as America, South America, Europe, Southeast Asia, The Middle East etc.