

**Rincian Produk:**

1. Material: Non woven polypropylene 100%
2. Size:80\*190+20cm,90\*200+20cm,160\*190+20cm,160\*200+20cm,180\*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; Putih ; other size available
5. Packing: 10pcs/bag, 100pcs/carton

## Indications□

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

## Specifications:

Bahan	Non woven polypropylene 100%
Warna	Putih biru
Gaya	elastics on width
Ukuran	80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
sedang mengemas	10 pcs/box, 100pcs/ctn
Bobot	25-65gsm or customized
Kemasan	PE bag coated, then paper carton package or customized
Fitur	Durable, soft, disposable
Penggunaan	SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

**Pelayanan kami:**

1. Provide the highest quality products  
Professional QC team strictly control each step of the production until loading.
2. The fastest delivery time  
If urgent, 1\*40HQ can be finished within 5-7days!  
Normally, 1\*40HQ can be finished within 10-15 days!  
Busy season, 1\*40HQ can be finished within 25-30 days!
3. The most reasonable price  
Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!
- 4.The sweet and quick service  
Our workers will be always so patiently making a response of your any questions within 5-10 minutes.
- 5.We can also help to source connected products.  
Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

