

Disposable gown is made from high-performance isolation protective material.

Widely used in medical safety protective, aseptic workshop, protective isolation, mining electronics factory, food factory farm animal husbandry biohazard and so on.

With two kinds processes: sewing and ultrasound thermal bonding.

Material features :

1. PP disposable gown

It is made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Color: white, green, blue, yellow, orange and so on

Material weight: 20-65gsm.

2. PP+PE disposable gown

It is made from PP+PE material, Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

Color: white, green, blue, yellow.

Material weight: 40-65gsm.

3. SMS disposable gown (Standard and Reinforced):

It is made from hydrophobic SMS/PP/PP+PE material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: white, blue.

Material weight: 35-65gsm.

Disposable gown style:

Tie/ hook and loop fastener on collar, tie on waist, elastic/knitted cuff

Disposable gown Specifications

| | |
|-------------|---|
| Description | disposable level1 level2 level3 isolation gowns SMS/PP+PE/PP protective cloth |
| Material | PP, PP+PE, SMS, CPE |

| | |
|----------------|--|
| Weight | PP 20-65gsm; PP+PE 40-65 gsm; SMS 35-65 gsm (comfortable,without glass fibres,latex free.) |
| Style | tie/ hook and loop fastener on collar,tie on waist,elastic/knitted cuff |
| Size | 120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm |
| Color | white,green,blue,yellow |
| packing | non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn |
| packing design | all inner box and carton printing can be as your design |
| carton size | 58X28X48 cm or as your requirement |
| Applications | hospital clinical medical personnel and patients. dust-free workshop,laboratory, food industry, electronic manufactures and so on |
| MOQ | 10000 pcs |
| Sample | can be supplied free for your quality checking within 3 days |
| Payment term | T/T,PayPal |
| Delivery time | as your quantity(we can finish 20,000 pcs per week) |
| FOB port | Wuhan/Shanghai/Ningbo/Guangzhou |



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

