

Beard cover surgery , high quality beard cover food service , Polypropylene beard nets

Specifications of beard cover surgery

Material	PP non woven
Fabric weight	Could be as your request
Ear loop	Elastic ear loop
Size	19Inch, 21 Inch, or as your requirements
Color	White,blue,green....etc
Packing	100 Pcs/PE bag ,1000 Pcs/carton, or as your requirement
Packing design	All inner box and outer carton printing can be as your design
Applications	Hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on
MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Packaging of beard cover food service

Specification	Normal packing	Carton size
19 Inch	100 Pcs/PE bag ,1000 Pcs/carton	26*22*32 cm
21 Inch	100 Pcs/PE bag ,1000 Pcs/carton	52*23*26 cm

Picture of Polypropylene beard nets



FAQ

Q1: .How do you solve quality issues if any?

A: Our factory will pay the highest attention to quality issues. We will deeply investigate the root cause of the problems and find the best solution to avoid any future problems. Meanwhile we will corresponding compensation for our customers.

Q2: What is your products expiration date?

A: Around 3-5 years.

Q3: Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.