

Indication of disposable dental gowns

- 1) Hospital, clinical, casualty wards, medical personnel and patients ;
- 2) Mining electronics factory, food factory, animal farm, husbandry bio-hazard industry to isolate dust...etc

We have following kinds of disposable dressing gowns

- 1) PP isolation gowns
- 2) PP+PE surgical gowns
- 3) SMS surgical gowns
- 4) CPE gowns

The Specification of disposable dental gowns

Material	PP+PE
Color	White, green, blue, yellow, etc.
Cuff	Knitted cuff or elastic cuff,
Neck	Ties or velcro
Waist	Tie on
Sleeves	Long sleeves/short sleeves/non sleeves
Sterilization	EO Sterilization or non-sterilize
Reinforced layer	as your requirements
Size	S-6XL, or as your request
Minimum order quantity	5,000 pcs
Price:	Negotiable
Payment terms:	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability:	1*40 HQ can be finished in 15 days
Delivery time:	According to the quantity, normally, 30 days after confirmed the order

Packing of disposable dressing gowns

Sterile packing	1 pcs/ sterile paper plastic bag, 50 pcs/Ctn
Regular packing	5 pcs/bag, 50 pcs/Ctn
Remark: all can be special customized according to your actual requirements	

The size chart of the disposable exam gowns

Size	S	M	L	XL
Height (CM)	110	115	120	125
chest circumference(CM)	30	137	140	150

The picture of the disposable dressing gowns





FAQ:

(1) Are you a manufacturer or a trading company?

We are a manufacturer located in Xiatao city, Hubei. Welcome to visit our factory.

(2) May I get samples to test quality? Small trail order acceptable?

Yes, we'd like to send you samples for evaluation. Small trail order is also welcomed.

(3) Can you offer customized service?

Yes, customized size and OEM brand available.

(4) What's your main products?

Our main products include disposable face mask, disposable bouffant cap, disposable coverall, lab coat, surgical gown, disposable cap, apron, shoe cover, etc.

(5) What's your delivery time?

Usually within 15 days for 20feet container.