

Property Disposable PP lab suit

All PP dress. Made of high quality material, providing excellent protection and superior comfort:

1. Size: S-5XL Or customized
2. Made of Nonwoven Material, Breathe And comfortable, clean and sanitary.
3. Comfortable and light, easy to do.

Specifications Disposable PP lab suit

Name.	Disposable PP lab coat.
Material	PP.
Style.	Elastic cuff or knitted
Size	S-5XL
Color	Blue
Packing.	10. Pcs / Box, 10 Bags / CTN
Packing design	All boxes in and carton printing can be your design
Moq	5000 pcs
Samples	Can be provided free for your quality inspection inside 2. Days
Payment Method	T / T, Paypal, Western Union, L / C, D / A, D / P, Escrow,
Delivery time	Capacity: 10,000 pcs / day, 1 * 40HQ can be finished in 20. -25. Day;
Port FOB	Wuhan / Shanghai / Ningbo / Guangzhou

Picture Disposable PP lab coat.







