

All [disposable reinforced surgical gown](#) are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

|                        |  |
|------------------------|--|
| <b>Fabric Material</b> | SMS material   |
| <b>Style</b>           | With Velcro on collar , tie on waist knitted cuffs             |
| <b>Dimensions size</b> | X-3XL  |
| <b>Fabric weight</b>   | 30-60g, could be as your request                               |
| <b>Color</b>           | Blue color   |
| <b>Certificate</b>     | ISO13485, FDA  |
| <b>MOQ:</b>            | 10000 pcs  |
| <b>Packing design</b>  | can be as your design  |
| <b>Packing</b>         | 1pc/plastic paper bag,50pcs/ctn, or as your requirement        |
| <b>Capability</b>      | 60,000pcs per day  |
| <b>Delivery time</b>   | 7-30 days  |
| <b>Applications</b>    | laboratory, clinic, food industry, dust-free workshop,hospital |

The picture of [disposable reinforced surgical gown](#)



Set-in Sleeves



Raglan Sleeves



Waist card, more  
nient distinction the belt



Heat-welded ties



Ultra-sonic heat  
sealing seams





