

Surgical PP Nonwoven Disposable Round Bouffant Caps Manufacturers

Specification:

Item	Surgical PP Nonwoven Disposable Round Bouffant Caps Manufacturers
Material	9g, 10g, 12g, 14g, 18g, 20g, 25g etc
Style	Elastic (latex free), single and double elastic style
Color	white, blue, green
packing	100pcs/bag, 1000pcs/ctn
packing design	All inner box and carton printing can be as your design
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS
Quality Control	1.We will not start to make the products until you confirm the sample. 2.The packing will be water-proof, damp-proof and sealed
Certificate	CE, ISO13485, FDA





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL



