

disposable nonwoven surgical caps for doctor

Specification:

Material	PP, SPP, SBPP, SMS, SMMS, Microporous
Color	White, blue, green, yellow, red....etc
Fabric weight	PP 10-30gsm, SMS 30-65gsm, Microporous45-65gsm
Type	ties on
Elastic	Single/double
Size	18",19",20",21",24"....etc
Certificate	CE,ISO13485,FDA
Packaging	100/bag,2000pcs/ctn, or as your requirements
Sample	Sample can be finished in 2 days
Model	disposable surgical cap

Advantages:

1. Directly manufacturer, disposable mop cap, any colors.
2. ISO, CE, good protect, latex free, easy breath, eco friendly,
3. Capacity: 200,000pcs/day, 1*40HQ can be finished within 10-15 days;
4. OEM packing design, service in 24 hours.