

## Product details

### Introduction of plastic blue CPE isolation gown:

CPE disposable isolation gown is made from hydrophobic CPE material and can be used in hospital, laboratory, food industry, or any where need protection. Weight can be made from 30-65g. And CPE material with the features of Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

### Specification of plastic blue CPE isolation gown:

material	CPE
Weight	CPE 30-65 gsm; (comfortable, without glass fibres, latex free.)
Style	Neck: all in one Waist: ties on Cuff: with thumb loop Opening on back
Size	120x130cm, 130x150cm, 150x160cm, 150x190cm, 170x170cm
Color	blue, or yellow, pink, red, grey, black etc.
packing	non sterile, 10 pcs/bag, 10 bags/ctn. sterile, 1 pcs/sterile bag, 100 pcs/ctn
packing design	All the inner box and carton printing can be as your design
carton size	Standard packing carton 58X28X48 cm
Applications	hospital clinical medical personnel and patients. used in hospital, laboratory, food industry, and any where need protection

## Product display




orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

- 1, Free sample.
- 2, Small order acceptable;
- 3, L/C, T/T, D/P terms available;
4. FOB, CNF, CIF, EXW terms available;
- 5, Competitive and stable price because we are factory.
- 6, Qualified products. 3rd part product inspection is acceptable.
- 7, OEM service. More options on material, pattern, size and packing;
- 8, Timely delivery. The best lead time is 25-30 day for a 40HQ container.
- 9, Professional export service. Documentation, clearance, sea transportation, all in one.
- 10, Good after sale service. We take care complaints as we would like to be your long-term supplier!