

**Product Details:**

1. Material: Non woven with PE filme(microporous)
2. Size:80\*190+20cm or customized size
3. Elastic: 5mm along side width+30cm on lenght
4. Colour: Blue; White ; other size availabe
5. Packing: 10pcs/bag, 100pcs/carton

**Indications**□

- 1, Widely used for hospital,SPA, Hotel, Massage, Salon, one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

## Specifications:

Material	Non woven with PE filme(microporous)
Color	White, blue
Style	elastics on width
Size	80*190+20cm,or customized size
packing	10 pcs/box, 100pcs/ctn
Weight	70gsm or customized
Package	PE bag, then carton package or customized
Feature	Durable, soft, disposable
Usage	Hospital, SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

**Our Services:**

1. Provide the highest quality products

Professional QC team strictly control each step of the production until loading.

2. The fastest delivery time

If urgent, 1\*40HQ can be finished within 5-7days!

Normally, 1\*40HQ can be finished within 10-15 days!

Busy season, 1\*40HQ can be finished within 25-30 days!

### 3. The most reasonable price

Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!

### 4.The sweet and quick service

Our workers will be always so patiently making a response of your any questions within 5-10 minutes.

### 5.We can also help to source connected products.

Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!

<https://www.china-disposable.com/products/disposable-fitted-sheet.html#.YLnLuvkzaUk>