

La proprietà di Disposable SMS lab coat

Tutti cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: S-5XL or customized
2. Made of nonwoven Materiale, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

The Specification of Monouso sms Camice da laboratorio

Nome	Monouso sms Camice da laboratorio
Materiale	sms
Stile	Knitted collar and knit cuffs, with bottons and pockets
Misurare	S-5XL
Colore	blu
Imballaggio	5pcs/box,10 bags/ctn
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Campione	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western Union, L/C,D/A,D/P,ESCROW,
Tempo di consegna	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

La proprietà di Disposable SMS lab coat

Tutti cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: S-5XL or customized
2. Made of nonwoven Materiale, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

The Specification of Monouso sms Camice da laboratorio

Nome	Monouso sms Camice da laboratorio
Materiale	sms
Stile	Knitted collar and knit cuffs, with bottons and pockets
Misurare	S-5XL
Colore	blu
Imballaggio	5pcs/box,10 bags/ctn
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Campione	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western Union, L/C,D/A,D/P,ESCROW,
Tempo di consegna	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

