

Dettagli del prodotto:

1. Material: Non woven polypropylene 100%
2. Size:80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; Bianco ; other size available
5. Packing: 10pcs/bag, 100pcs/carton

Indications□

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

Specifiche:

Materiale	Non woven polypropylene 100%
Colore	Bianco blu
Stile	elastics on width
Misurare	80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
Imballaggio	10 pcs/box, 100pcs/ctn
Il peso	25-65gsm or customized
Pacchetto	PE bag coated, then paper carton package or customized
Caratteristica	Durable, soft, disposable
Utilizzo	SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

I nostri servizi:

1. Provide the highest quality products
Professional QC team strictly control each step of the production until loading.
2. The fastest delivery time
If urgent, 1*40HQ can be finished within 5-7days!
Normally, 1*40HQ can be finished within 10-15 days!
Busy season, 1*40HQ can be finished within 25-30 days!
3. The most reasonable price
Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!
- 4.The sweet and quick service
Our workers will be always so patiently making a response of your any questions within 5-10 minutes.
- 5.We can also help to source connected products.
Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

