

XXXXXXXXXXXXXXXXXXXXXXXX

尺寸	规格	数量	1 * 40 HQ 箱装数量 PCS
17.5 * 9.5cm 3-层 18 + 20 + 20g	50 片 / 40 箱 / CTN	52 * 19 * 42cm	3,280,000 片 PCS
17.5 * 9.5cm 3-Ply 25 + 25 + 25g	50 片 / 40 箱 / CTN	51 * 37 * 30cm	2,400,000 片

XX

XXXXXXXXXXXXXXXXXXXXXXXX

品牌	Xiantao 湖北
材质	PP Non Woven + Filter + PP
规格	15 + 20 + 25g / 20 + 20 + 25g / 20 + 25 + 25g
包装	每箱 40 箱
尺寸	17.5x9.5cm / 14.5x9.5cm / 12x7cm
用途	医用口罩
重量	每箱重量
外包装尺寸	52x38x30cm / 52x19x42cm
内包装尺寸	每包尺寸

XXXXXXXXXXXXXXXXXXXXXXXX







 orient Honest  
MEDICAL INTERNATIONAL

