

4ply KF94

- linting
-
-

4ply KF94

material	1 PP Nonwoven 2 3 4 PP
	4ply
s ize	20*8 cm
w 8	30 GSM + 25GSM + 25GSM+25GSMc be
color	
	ce iso13485 fda
MOQ	20000pcs
	5 0 pcs/b ox 200 0 PC s/ctn
	100 1PC
	7-30

4ply KF94







orient Honest
MEDICAL INTERNATIONAL



orient Honest
MEDICAL INTERNATIONAL

