



Model	3D KF94
BFE/PFE	BFE ≥ 95%
Material	PP + Nonwoven + Nonwoven + PP
Color	White and Brown
Weight	50g + 25g + 25g + 30g
Size	20x8cm
Quantity	50 pcs/bag, 40 bags/carton, 2000 pcs/lot
Lead Time	7-15 days
Payment	T/T, L/C, X-Transfer, PayPal
Service	OEM/ODM available
Standard	GB2626-2019

orient MEDICAL INTL

✔ With high moisturizing transferability technique, it absorbs moisture quickly, boasting a fresh fit even if worn for a long time, and prevents misting even when wearing glasses.



* 4層構造

* 3D構造

* 高保湿転送技術

* 快適

* 長時間着用可能

マスク 高保湿転送技術











