

Product details

1. Material: PP

2. Color: White

3. Capacity: 10,000 pieces

4. Weight: 1.5kg

5. Size: 10cm x 10cm x 10cm

Material: **PP**

CPE: 100%

PP / PE / SMS / SMMS / PP

Weight: 1.5kg

Size: S-6XL

Capacity: 10,000 pieces

Weight: 10kg / CTN

Weight: 1.5kg

MOQ: 5000

Weight: 3kg

T / T / L / C / D / A / P

Capacity: 10,000 / 1 * 40HQ 10 15

FOB: 1.5 / 1.5 / 1.5

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a factory?

A: Yes, we are a factory established in 2003. We have a professional R&D team and a mature production line.

Q2. What is your price range?

A: Our price range is from USD 18 to 20 per unit.

Q3. Do you have any certifications?

A: Yes, we have ISO 9001, CE, ISO 13485, FDA, REACH, and NELSON certifications.

Q4. Can you provide OEM services?

A: Yes, we can provide OEM services with a lead time of 3-7 days.