

Product details

2

- PP2&038;V
- 1.2&038;
- 2.5
- 3.No
- 4.PP;
- 5.
- 6.
- 7.

2

	S-5XL
	PP50g
	50
MOQ	5,000PCS
	1. With2&038; 2. Elastic5 3. No
	1suit / 50suits / CTN
	1 * 40HQ 15-20days
	T / TL / CD / AD / P
	2
	ISO13485 CEFDAType5&038; 6

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a manufacturer?

A: Yes, we are a manufacturer since 2003. We have a factory in China with a production capacity of 2013 units per month.

Q2. What is your minimum order quantity (MOQ)?

A: Our MOQ is 100 units, with a price range of USD 18 to 20 per unit.

Q3. Do you have any certifications?

A: Yes, we are certified with ISO 9001, EN 14001, CE, and REACH. We also have FDA and NELSON certifications.

Q4. Can you provide samples?

A: Yes, we can provide samples. The lead time for samples is 3-7 days.

The happy time with customers

