

## Product details

1. Material: Spunbonded

2. Weight: 25gsm

3. Size: S - 6XL

4. Quantity: 10 pieces / 100 pieces / CTN

5. Usage:

6. Usage: For home use, hospital use, etc.

7. Usage: For home use, hospital use, etc.

8. Usage: For home use, hospital use, etc.

Material	PP
Weight	PP25g
Size	S-6XL
Quantity	10 pieces / 100 pieces / CTN
Usage	For home use, hospital use, etc.
Usage	For home use, hospital use, etc.

## Product display




orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q1. Are you a factory?**

A: Yes, we are a factory since 2003. We have a professional R&D team and a production line since 2013. We can provide you with the most competitive price.

**Q2. What is your price range?**

A: Our price range is USD 18-20 per unit.

**Q3. Do you have any certifications?**

A: Yes, we have CE, ISO 13485, FDA, REACH, and NELSON certifications.

**Q4. Can you provide a sample?**

A: Yes, we can provide a sample within 3-7 days.