

## Product details

Product Name: [Placeholder]

1. Material: PP50g
2. Color: [Placeholder]
3. Weight: 50gsm
4. Thickness: 2mm / [Placeholder]
5. Size: S-2XL
6. Quantity: [Placeholder]

Product Description: [Placeholder]

Material	PP50g [Placeholder]
Weight	PP50g
Thickness	2mm [Placeholder]
Size	S M L XL 2XL 3XL [Placeholder]
Color	[Placeholder]
Quantity	1-5pcs / [Placeholder] 10 [Placeholder]
Product Description	[Placeholder]
MOQ	5000 [Placeholder]
Lead Time	3 [Placeholder]
Payment Terms	T / T [Placeholder] L / C [Placeholder] / A [Placeholder] / P [Placeholder]
Price	[Placeholder] 10,000 [Placeholder] 1 40HQ [Placeholder] 10 [Placeholder] 15 [Placeholder] *
FOB [Placeholder]	[Placeholder] / [Placeholder] / [Placeholder]

## Product display



orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q1. Are you a manufacturer?**

A: Yes, we are a manufacturer since 2003. We have a factory in China with a production capacity of 2013 units per month.

**Q2. What is your price range?**

A: Our price range is from USD 18 to 20 per unit.

**Q3. Do you have any certifications?**

A: Yes, we have ISO 9001, CE, ISO 13485, FDA, REACH, and NELSON certifications.

**Q4. Can you provide samples?**

A: Yes, we can provide samples. The lead time for samples is 3-7 days.

