

Product details

1. Material

2. Made in

3. Grade

4. Widely used in

1. Material

2. Made in

3. Grade

4. Widely used in

5. Widely used in

6. Widely used in

1. Material

Material	disposable
Material	LDPE / CPE
Thickness	5mm - 25mm
Weight	0.01 - 0.045g
Color	
Size	* 107mm x 61mm x 71mm * 117mm x 70mm x 80mm * 140mm x 70mm x 80mm
Quantity	
MOQ	1000 / 10 boxes / CTN
Lead time	
MOQ	10 PCS
Lead time	3 days
Lead time	T / T / L / C / D / A / D / P /
Lead time	1 * 40HQ / 25
FOB	

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a factory?

A: Yes, we are a factory. We have been established since 2003 and have a production capacity of 2013 units per month.

Q2. What is your price range?

A: Our price range is USD 18 to 20 per unit.

Q3. Do you have ISO 9001 and CE certification?

A: Yes, we have CE, ISO 13485, FDA, REACH, and NELSON certifications.

Q4. Can you provide samples?

A: Yes, we can provide samples. The lead time for samples is 3-7 days.

