

Product details

Product description

3. Handmade kintted

2. Non-stretchable fabric

Product features

5 PC / 10bags / 1/50/

Product specifications

Size chart

Unit: cm

Product specifications

CM	S	M	L	XL	XXL	XXXL	XXXXL	XXXXXL
Length	100	103	107	110	112	113	120	120
Length	132	135	137	140	142	147	150	158

Product details

Material: coatm CE ISO

Material: PP

Weight: 50gsm

Product features: / / / /

Product features: workshop

Product features:

MOQ: 5000

Product features: 50/100/ CTN

Product features: 80

Product features: T / T L / C D / A D / P

Product features: OEM

Product features: CE ISO13485

Product features: Capacity 10,000 1 * 40HQ 10 15

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a manufacturer or a trading company?

A: We are a manufacturer since 2003. We have our own factory and R&D department.

Q2. What is your minimum order quantity (MOQ)?

A: Our MOQ is 1000 units. For some products, we can accept orders as low as USD 1800-2000.

Q3. Do you have any certifications or standards?

A: Yes, we are certified by CE, ISO 13485, FDA, REACH, and NELSON.

Q4. Can you provide samples for testing?

A: Yes, we can provide samples for testing. The cost of samples is 3-7 USD per unit.