

1. 品牌: ORIENT HONEST
 2. 型号: OH-1001 / OH1002 / OH1003
 3. 认证: CE / ISO13485 / FDA
 4. 材质: PP
 5. 重量: 10g-30g
 6. 尺寸: 18" x 19" x 21"
 7. 包装: 100 PC / 20 PC/CTN

品牌	ORIENT HONEST
型号	OH-1001 / OH1002 / OH1003
认证	CE / ISO13485 / FDA
材质	PP
重量	10g-30g
尺寸	18" x 19" x 21"
包装	100 PC / 20 PC/CTN
品牌	ORIENT HONEST
型号	OH-1001 / OH1002 / OH1003
认证	CE / ISO13485 / FDA
材质	PP
重量	10g-30g
尺寸	18" x 19" x 21"
包装	100 PC / 20 PC/CTN
品牌	ORIENT HONEST
型号	OH-1001 / OH1002 / OH1003
认证	CE / ISO13485 / FDA
材质	PP
重量	10g-30g
尺寸	18" x 19" x 21"
包装	100 PC / 20 PC/CTN
品牌	ORIENT HONEST
型号	OH-1001 / OH1002 / OH1003
认证	CE / ISO13485 / FDA
材质	PP
重量	10g-30g
尺寸	18" x 19" x 21"
包装	100 PC / 20 PC/CTN

1. 品牌: ORIENT HONEST

1. 品牌: ORIENT HONEST
2. 型号: OH-1001 / OH1002 / OH1003
3. 认证: CE / ISO13485 / FDA
4. 材质: PP
5. 重量: 10g-30g
6. 尺寸: 18" x 19" x 21"

□□□□□□□□□□□□□□□□









