

### The property of Disposable White PP lab coat

Our cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

1. Size: S-5XL or customized
2. Made of nonwoven fabric, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight, easy to carry out.

### The Specification of Disposable White PP lab coat

|                |                                                                    |
|----------------|--------------------------------------------------------------------|
| Product        | Disposable White PP lab coat                                       |
| Material       | pp                                                                 |
| Feature        | Knitted collar and knit cuffs, with buttons and pockets            |
| Size           | S-5XL                                                              |
| Color          | White                                                              |
| Quantity       | 5pcs/box, 10 bags/ctn                                              |
| Packing design | all inner box and carton printing can be as your design            |
| Moq            | 5000 pcs                                                           |
| Lead time      | can be supplied free for your quality checking within 2 days       |
| Payment way    | T/T, PayPal, Western Union L/C, D/A, D/P, ESCROW,                  |
| Capacity       | Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days; |
| FOB port       | Wuhan/Shanghai/Ningbo/Guangzhou                                    |

### The picture of Disposable White PP lab coat





