

□□□□

- 1. Material: Non woven polypropylene 100%
- 2. Size:80\*190+20cm,90\*200+20cm,160\*190+20cm,160\*200+20cm,180\*200+20cm or customized size
- 3. Elastic: 5mm along side width+30cm on lenght or customized size
- 4. Colour: Blue;□ ; other size availabe
- 5. Packing: 10pcs/bag, 100pcs/carton

Indications□

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

□□

|       |                                                                                   |
|-------|-----------------------------------------------------------------------------------|
| □□    | Non woven polypropylene 100%                                                      |
| □     | □□                                                                                |
| □□□□  | elastics on width                                                                 |
| □□□   | 80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size |
| □□    | 10 pcs/box, 100pcs/ctn                                                            |
| □□    | 25-65gsm or customized                                                            |
| □□□□□ | PE bag coated, then paper carton package or customized                            |
| □□    | Durable, soft, disposable                                                         |
| □□□   | SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home                             |

□□□□□□□□

- 1. Provide the highest quality products  
Professional QC team strictly control each step of the production until loading.
- 2. The fastest delivery time  
If urgent, 1\*40HQ can be finished within 5-7days!  
Normally, 1\*40HQ can be finished within 10-15 days!  
Busy season, 1\*40HQ can be finished within 25-30 days!
- 3. The most reasonable price  
Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!
- 4.The sweet and quick service  
Our workers will be always so patiently making a response of your any questions within 5-10 minutes.
- 5.We can also help to source connected products.  
Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

