



Model	3D KF94
BFE/PFE	BFE≥95%
Material	PP+Nonwoven+Nonwoven+pp
Color	Green/White
Weight	50g+25g+25g+30g
Size	20x8cm
Quantity	50pcs/bag, 40bags/carton, 2000pcs/lot
Lead time	7-15 days
Payment	T/T, L/C, X-Transfer, PayPal
Service	OEM/ODM
Standard	GB2626-2019

orient HOME MEDICAL INTERNATIONAL

✔ With high moisturizing transferability technique, it absorbs moisture quickly, boasting a fresh fit even if worn for a long time, and prevents misting even when wearing glasses.



* 4-layer construction

* 3D micro-ventilated design for breathability

* Soft and comfortable fit

* Effective filtration

* Durable and reusable

Product Name: [Placeholder]









