



Model	3D KF94
BFE/PFE	BFE≥95%
Material	PP+Nonwoven+Nonwoven+ pp
Color	Black, Green, White
Weight	50g+25g+25g+30g
Size	20x8cm
Quantity	50pcs/bag, 40bags/carton, 2000pcs/lot
Lead time	7-15 days
Payment	T/T, L/C, X-Transfer, PayPal
Service	OEM/ODM
Standard	GB2626-2019

orient HO MEDICAL INTERIORS

✔ With high moisturizing transferability technique, it absorbs moisture quickly, boasting a fresh fit even if worn for a long time, and prevents misting even when wearing glasses.



- * 4-layer construction
 - * 3D micro-ventilated design for breathability
 - * Soft and comfortable fit
 - * Effective moisture management
 - * Durable and long-lasting
- Additional features and benefits:











