

[korean KF94 4 Ply for Adults kf94](#)

Specifications

- linting
-
-

[korean KF94 4 Ply for Adults kf94](#)

material	1 PP Nonwoven 2 3 4 PP
	4ply
s ize	20*8 cm
w 8	30GSM + 25GSM + 25GSM+25GSM c be
color	
	ce iso13485 fda
MOQ	20000pcs
	5 0 pcs/b ox 200 0 PC s/ctn
	100 1 PC
	7-30

[korean KF94 4 Ply for Adults kf94](#)



orient Honest
MEDICAL INTERNATIONAL



orient Honest
MEDICAL INTERNATIONAL





 orient Honest
MEDICAL INTERNATIONAL

