

The property of Disposable SMS lab coat

Our cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: S-5XL or customized
2. Made of nonwoven fabric breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

The Specification of Disposable SMS

Item	Disposable SMS
Material	SMS
Collar	Knitted collar and knit cuffs, with buttons and pockets
Size	S-5XL
Weight	
Packing	5pcs/box,10 bags/ctn
Packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Lead time	can be supplied free for your quality checking within 2 weeks
Payment way	T/T, PayPal, Western Union L/C,D/A,D/P,ESCROW,
Capacity	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The property of Disposable SMS lab coat

Our cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: S-5XL or customized
2. Made of nonwoven fabric breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

The Specification of Disposable SMS

Item	Disposable SMS
Material	SMS
Collar	Knitted collar and knit cuffs, with buttons and pockets
Size	S-5XL
Weight	
Packing	5pcs/box,10 bags/ctn
Packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Lead time	can be supplied free for your quality checking within 2 weeks
Payment way	T/T, PayPal, Western Union L/C,D/A,D/P,ESCROW,
Capacity	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

