

Product details

Specifications

- 1. SMS type: 2D Barcode
 - 2. Weight: 20~100gsm
 - 3. Size: 100mm x 150mm
 - 4. Material: 100% Recycled Paper
 - 5. Dimensions: 200cm * 120cm / 200cm * 90cm / 180cm * 120cm / 180cm * 90cm / 90 * 60cm
- Minimum order quantity: 1000000

Usage

- 1. For SMS verification
- 2. For SMS authentication

SMS

Quantity	1 / 50 PCS / CTN
Weight	5 / 50 PCS / CTN
Minimum order quantity: 1000000	
Material: 100% Recycled Paper	
Size: 100mm x 150mm	
Weight: 20~100gsm	
Dimensions: 200cm * 120cm / 200cm * 90cm / 180cm * 120cm / 180cm * 90cm / 90 * 60cm	

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

