

1. 18mm x 21mm x 24mm

2.

3.

4.

5.

1.

2. FDA CE ISO13485

3.

4.

5.

6.

7.

8.

9.

Item	Quantity	Dimensions (CM)	N.G (kg / ctn)	W.G (kg / ctn)
100mm / 2000pcs / ctn	200,000	36X26X48	7.5	8
	1,150,000 1X20GP	36X26X48	7.5	8
	2,400,000 1X40GP	36X26X48	7.5	8
	2,800,000 1X40HQ	36X26X48	7.5	8

Why choose us?



SELLING POINT



SELLING POINT



SELLING POINT




SELLING POINT



01

Elastic design High Quality



02

High quality non woven cap



QUALITATIVE ANALYTICAL

