

The feature of surgical lab coat

- ü Isolate external liquids
- ü Filter fine particles and dust
- ü Comfortable and breathable
- ü Suitable for men and women

ü Environmental protection fabric

Fabric Material	SMS
Style	knitted cuffs and collar ,with buttons
Dimensions size	S-3XL
Fabric weight	35-65g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	10-30 days
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

The picture of Surgical Lab coat with knitted cuffs medical supplier




orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

