

Product details

SMS
 PPE
 EC 686/89 / EEC
 3
 SMS
 5
 6
 686/89 / EEC

1
 2 CE
 CE 3
 6

SMS features
 1
 1
 2
 3 2
 4 Adhesive
 5
 6 OEM OPTION;
 7 CATEGORY iii
 1 TYPE 5 EN 13982 1 2004 airborne PPE
 2 TYPE 6 EN 13034 + EN 468 PPE
 SMS 8

SMS

	SMS 2 SMS
	S-5XL
	SMS45g
	45
MOQ	5,000PCS
	BLUE
	1. With 2 & 038; 2. Elastic 5 3. No
	1suit / 50suits / CTN
	1 * 40HQ 15-20days
	T / TL / CD / AD / P
	2
	ISO13485 CE FDA Type5 & 038; 6

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a manufacturer or a trading company?

A: We are a manufacturer since 2003. We have our own factory and R&D department.

Q2. What is your minimum order quantity (MOQ)?

A: Our MOQ is 1000 units. For some products, we can accept orders as low as USD 1800-2000.

Q3. Do you have any certifications or standards?

A: Yes, we are certified with CE, ISO 13485, FDA, REACH, and NELSON.

Q4. Can you provide samples for testing?

A: Yes, we can provide samples. The lead time for samples is 3-7 days.

The happy time with customers

