

MOQ 5,000pcs

MOQ	5,000pcs
Material	PE
Size	90 * 120mm / 102 * 128mm
Weight	0.015g
Color	White
Surface	Smooth
Thickness	0.1mm
Quantity	1000 / OPP / 250pcs / CTN / 51 * 30 * 23cm / CTN
Lead time	11 days / CTN

MOQ 5,000pcs

MOQ 5,000pcs

MOQ 5,000pcs

MOQ 5,000pcs

1. Our PE material is made of high quality polyethylene material. The thickness is 0.1mm. The weight is 0.015g. The size is 90 \* 120mm / 102 \* 128mm. The color is white. The surface is smooth. The quantity is 1000 / OPP / 250pcs / CTN / 51 \* 30 \* 23cm / CTN. The lead time is 11 days / CTN.

2. The material is made of high quality polyethylene material. The thickness is 0.1mm. The weight is 0.015g. The size is 90 \* 120mm / 102 \* 128mm. The color is white. The surface is smooth. The quantity is 1000 / OPP / 250pcs / CTN / 51 \* 30 \* 23cm / CTN. The lead time is 11 days / CTN.

MOQ 5,000pcs



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

