

Product details

1. SMS 50g

1. Softable ability
2. Prevent TE
3. Prevent TE
4. Prevent TE

2. SMS 50g

1. Prevent TE
2. Prevent TE
3. Prevent TE

3. SMS 50g

1	SMS50g
2	Prevent TE
3	120x130cm 130x150cm 150x160cm 150x190cm 170x170cm
4	Prevent TE
5	58X28X48
6	10/10/ CTN 1/100/ CTN
7	Prevent TE
8	PP 15-65gsm PP + PE 40-65 GSM SMS 35-65 GSM
9	Prevent TE

4. SMS 50g

1. Prevent TE
2. ISO13485 CE
3. 10,000/1 * 40HQ 10-15days finisehd
4. OEM 24
5. 3

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a manufacturer or a trading company?

A: We are a manufacturer since 2003. We have a factory in China with a production capacity of 2013 units per month.

Q2. What is your minimum order quantity (MOQ)?

A: Our MOQ is 100 units. The price is USD18 per unit, with a 20% discount for orders over 1000 units.

Q3. Do you have any ISO or CE certification?

A: We have CE, ISO13485, FDA, REACH, and NELSON certifications.

Q4. Can you provide a sample?

A: Yes, we can provide a sample. The sample fee is 3-7 units, which can be refunded upon bulk order.

The happy time with customers

