

**KF94**

- linting
- 
- 

**KF94**

<b>material</b>	1 PP Nonwoven 2 3 4 PP
	4ply
<b>s ize</b>	20*8 cm
<b>w 8</b>	30GSM + 25GSM + 25GSM+25GSM c be
<b>color</b>	
	ce iso13485 fda
<b>MOQ</b>	20000pcs
	5 0 pcs/b ox 200 0 PC s/ctn
	100 1 PC
	7-30

**KF94**



orient Honest  
MEDICAL INTERNATIONAL



orient Honest  
MEDICAL INTERNATIONAL





 orient Honest  
MEDICAL INTERNATIONAL

