

1. 100% cotton
2. 11*62 cm
3. 100pcs / bag
- 4.
5. 10-30 GSM
- 6.

Color

- 1.
2. 1
- 3.
- 4.

100pcs / bag 1000 pcs / ctn

Color

1. OEM
2. OEM
- 3.
- 4.

Port

Color	
Color	
Weight	10-30 GSM
Material	100% cotton
Size	11*62 cm
Material	
Quantity	100pcs / bag 1000 pcs / ctn
Moq	100000 PC
Quantity	100pcs / bag 1000 pcs / ctn
Quantity	3
FOB Port	Wuhan/Shanghai/Ningbo/Guangzhou