

Anti-bacteria 100% effective, 3-layer non-woven fabric
 100% cotton fabric

100% cotton fabric

- 1) 100% cotton fabric
- 2) 100% cotton fabric, size: 17.5 cm x 9.5 cm X
- 3) 100% cotton fabric, 100% cotton fabric
- 4) 100% cotton fabric / 100% cotton fabric / 100% cotton fabric / 100% cotton fabric / 100% cotton fabric / 100% cotton fabric, 100% cotton fabric 100% cotton fabric
- 5) 100% cotton fabric, ISO, CE, BFE 99%, 100% cotton fabric / 100% cotton fabric
- 6) 100% cotton fabric, 100% cotton fabric, 100% cotton fabric
- 7) 100% cotton fabric: 100% cotton fabric, 100% cotton fabric, 100% cotton fabric, 100% cotton fabric

100% cotton fabric

50PCS / box, 40boxes / CTN, 52x38x30cm;
 1 * 20 GP container 1,000,000 PCS per container;
 1 * 40 HC container 2,450,000 PCS per container.

100% cotton fabric per container.

100% cotton fabric

MOQ	100
Material	100% cotton fabric
Color	T / T, L / C, D / A, D / P, 100% cotton fabric, 100% cotton fabric, 100% cotton fabric
Quantity	1,000,000 / container, 1 * 40 HQ container 15 boxes per container
Weight	100% cotton fabric per container
Lead time	100% cotton fabric, 2 weeks per container

100% cotton fabric



orient Honest
MEDICAL INTERNATIONAL

