

## Product details

CE / ISO13485 PP

Material: PP

Weight: 28gsm

Size: S, M, L, X- 2X - , - , 4X -

Color: White

Packaging: polybag 10, 5 polybags

1 polybag, 50 polybags

Usage: , ,

**CE / ISO13485 PP**

Material	CE / ISO13485 PP, PP28g
Material	PP
Color	White
Weight	PP28g
Size	S-6XL
Packaging	10PCS / polybag, 100PCS / CTN
Usage	
Notes	

Notes:

1. Lead time: 3-5 days.

2. Minimum order quantity: 1000 pieces.

3. Payment terms: T/T

4. Delivery terms: FOB.

□□ □□

1. □□□ □□□ □□□ □□□ □□□ □□ □□ □□□□ □□□□.

2. □□□ □□, □□ □□ □ □□□□

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q : 1.** 貴社はどのような会社ですか、どのような製品を扱っていますか、教えてください。

A : 弊社は、主に中国産の電子部品を扱っており、2003年に設立されました。現在は2013年に比べて、より多岐にわたる製品を取り扱っています。

**Q : 2.** What is your minimum order quantity?

A : 通常は、ご注文金額が USD18 から 20 程度からとなります。

**Q : 3.** Do you have ISO, EN, CE, FCC, etc. certifications?

A : はい、ISO, EN, CE, FCC, ISO13485 などの認証を取得しています。REACH にも対応しています。

**Q : 4.** Can I get a sample? How much does it cost?

A : はい、サンプルは 3~7 日程度で提供可能です。送料は別途お見積りいたします。