

Product details

□ □□ □□ □□□ □□

- 1) □□: SMS, □□ □ □, □□ □□ □ □□ □□□ □□
- 2) □□ □□: 20-100gsm
- 3) □□: □□ □□
- 4) □□□: □□□ 2 □□ □, □ □□□□ □□□ □□ □ □□□ □□□ □□
- 5) □□: 200 cm * 120 cm, 200 cm * 90 cm, 180 cm * 120 cm, 180 cm * 90 cm, 90 * 60 cm
(□□□□ □□□ □□ □□ □ □□□□)

□□□ □ □ □□□ □□□ □□ □□ □□□□

- 1) □□□ □□ □ □□□ □□□ □□□ □□ □□ □□□□. □□□□□.
- 2) □□, □ □□ □□ □□ □ □ □□.

SMS □ □□□ □□ □□

□ □ □□	1 pc V □ □ □□ □□ □□, 50 pcV/Ctn
□□ □□	5 pcV□□, 50 pcV/Ctn
□: □□ □ □□ □ □□ □□ □□□ □□ □□ □□□	

Product display



orient Honest
MEDICAL INTERNATIONAL

