

5 6 0 000 000 00 :

1. 00 : 00 000 00 60g
2. 000 : 000 00 00, 00 00, 00, 00 000 000 00 0000.
3. 00 : 0 00, 0000 00 00 00 0 00 0000 00 0 0000
4. 00 : S-6XL
5. 00 : 000.
6. 00 : 1PC / polybag, 25pcs / CTN, 5PCS / polybag, 25pcs / CTN
1PCS / polybag, 50PCS / CTN, 10PCS / polybag, 100PCS / CTN.

5 6 0 000 000 00 :

00	00 000 00 00 00 0 50 00 00 6 0 000 00, 000 00 00 00
00	S : 96 * 124cm M : 100 * 130cm L : 107 * 136cm XL : 110 * 140cm 000 : 114 * 144cm
00	00 000 00
00	SF60G
000	000 00 00, 00 00, 00, 00 000 000 00 000
MOQ	5,000PCS
0	00
00	000 00 00 00, 00 000, 00 00, 00 00 00, 00 00 00 00 00 00 00 00000.
00	1PC / polybag, 25pcs / CTN, 5PCS / polybag, 25pcs / CTN 1PCS / polybag, 50PCS / CTN, 10PCS / polybag, 100PCS / CTN
00 00	1 * 40HQ 00000 15~20 00 00 0 0 0000.
00	00 00 00, 000 00, 000, 000, 000, 00.
00 00	T / T, L / C, D / A, D / P, 000 000, 0000, 000
00	00 00 2 0 000 0000 000 000 00 0 0 0000
000	ISO13485. CE, FDA

000 00 :

1. 0 00 00 00 0 00 0000 0000 0000 0000 00.
- 000 0 000 0000 2. 00.
3. 0000 0 0000 0000.
4. 00 000 0 00, 00 00 0 ..
- 00 00 5. 00
7. 000 00000 00 0000 00 0 0 00000.

00 000 00 00 0000 00 :

00, 00, 00, laboratorial, 00 00 1.

2. 000 00 0000 0000 000 00000, 00, 00, 00, 00, 00 0 0000 000 00.

00 000 00 00 0000 00 :



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

